

SurgiTel®

Internal Sales Bulletin: Surgical Mini LED

Loupe Configuration Options (TTL or FLM)



Brightness:	About 100,000 Lux (average)
Spot Size:	4.0" Diameter (10cm) at 20" (50cm)
Color Temperature:	5,700k
Beam Color:	Cool only
LED Life (average):	Up to 10,000 hours or over 5 years
Lens Diameter:	19mm
Weight (light unit & optics):	14.6g
Battery:	Rechargeable Lithium Ion Battery 7.4v 4,888 mAh
Battery Life:	Up to 12 hours
AC Charger:	Input 100-240v, 50-60Hz, 0.5A Output: +8.4v, 1.5A

Note: Only use the SHL850 Battery Control Pack with the Surgical Mini LED.