

Novaclip®

Bulldogs – Atraumatic Spring Clips

This product is single use, disposable and latex-free.

- Low profile due to Novaclip's® smaller design.
- Easy insertion and removal thanks to the attachment pads – Easy to apply manually or with surgical instruments.
- Novaclip® provides occlusion of normal or atherosclerotic vessels that are less than 4mm in diameter and may be used in a variety of applications.

Novaclip® surgical Spring Clips are indicated for use in:

- Peripheral surgery
- Vascular surgery
- Cardiovascular surgery
- General surgery



Full size view



N-10112



N-10113



N-10108

Reference	Jaws length	Jaws design	Jaws opening	Jaws angulation	Pressure	Color code	Packaging (per box)
N-10112	17 mm	Soft/Soft	12 mm	Straight	1/2 force	Lilac	15 units
N-10113	16 mm	Soft/Soft	12 mm	Angled	1/2 force	Yellow	15 units
N-10108	24 mm	Soft/Soft	17 mm	Straight	1/2 force	Blue	15 units

Applicator for Novaclip® spring Clips

Reference	Overall length
VSA 185	19 cm

