



Laparoscopic catheter suitable for Hemostats like 3-D Matrix PuraBond®

- The applicator exists of outer cylinder and soft inner tube
- When pushed out, the tip of the inner tube articulates
- Articulating tip enabling accurate application of the product



Technical information of the Laparoscopic catheter

		Specification
External diameter	5 mm	Fits in regular trocar port
Effective length	250 mm	Sufficient effective reach for most procedures
Material	Stiff outer cylinder with soft inner tube	When pushed out, the soft tip of the inner tube articulates, enabling accurate application of the product
Product reference	634163	Producer: TOP Corp., Tokyo, Japan

distributed by

medicovation
 Medicovation GmbH & Co KG
 Am Wiesenbusch 2
 45966 Gladbeck / Germany
 Phone: +49 2043 944 192

Ordering details

Packing Unit = 10 catheters in 1 Box; Minimum purchase Volume 1 Box;

Please send PO to Medicovation

Email: bestellung@medicovation.de

eFax: +49 2045 41 33 882