

# Greyhound®

## Bulldogs - Adjustable Spring Clips

**This product is single use, disposable and latex-free.**

The Greyhound® Adjustable Spring allows for variable occlusive force to meet the user needs, allowing small vessel occlusion.



### Adjustable to three occlusive force settings:

Setting window confirms the selected occlusive force.

- 1 = 1/2 force
- 2 = 3/4 force
- 3 = Full force



- 6mm jaw opening; compatible with 6mm Spring Clip Applicator
- Lightweight
- Low profile, soft jaw surface
- May be used for anchoring suture

Reference	Jaws length	Jaws design	Jaws opening	Jaws angulation	Pressure	Color code	Packaging (per box)
<b>N-10157</b>	<b>18 mm</b>	Soft	Adjustable	Straight	Adjustable 3 forces	Grey	15 units



Applicator for Greyhound® Spring Clips	
Reference	Overall length
<b>VSA 1805</b>	<b>18 cm</b>